

Xpert® BCR-ABL Ultra p190

For use with GeneXpert® System with Touchscreen



Catalog Numbers

REF RBCRABLP190-10

303-0947 | Rev. A | 2023-08

RUO For Research Use Only. Not for use in diagnostic procedures.

Trademark, Patents and Copyright Statements

Cepheid®, the Cepheid logo, GeneXpert®, and Xpert® are trademarks of Cepheid, registered in the U.S. and other countries. All other trademarks are the property of their respective owners.


THE PURCHASE OF THIS PRODUCT CONVEYS TO THE BUYER THE NON-TRANSFERABLE RIGHT TO USE IT IN ACCORDANCE WITH THESE INSTRUCTIONS FOR USE. NO OTHER RIGHTS ARE CONVEYED EXPRESSLY, BY IMPLICATION OR BY ESTOPPEL. FURTHERMORE, NO RIGHTS FOR RESALE ARE CONFERRED WITH THE PURCHASE OF THIS PRODUCT.

© 2020–2023 Cepheid.

See [Revision History](#) for a description of changes.

Table of Contents

	Getting Started.....	5
	Product Information.....	5
	Proprietary Name.....	5
	Common or Usual Name	5
	Product Description, Summary, and Principle of Procedure	5
	Product Description	5
	Summary and Explanation.....	5
	Principle of the Procedure.....	6
	Reagents, Instruments, and Materials.....	7
	Reagents	7
	Materials Provided.....	7
	Materials Required but Not Provided	8
	Materials Available but Not Required.....	8
	Warnings and Precautions	8
	Chemical Hazards, Storage and Handling	10
	Chemical Hazards.....	10
	Storage and Handling	11
	Specimen Collection, Testing, and Results	12
	Specimen Collection.....	12
	Specimen Collection, Transport and Storage	12
	Procedure	12
	Before You Start	12
	Preparing the Sample.....	12
	Preparing the Cartridge.....	13
	Starting the Test: GeneXpert System with Touchscreen.....	14
	Viewing Results: GeneXpert System with Touchscreen	15
	Quality Control	15
	Results.....	16
	Quantitative Results.....	16
	Troubleshooting Guide	19
	Retests.....	20
	Limitations	21
	Limitations of the Assay.....	21
	Specific Performance Characteristics	23
	Analytical Data	23
	Assay Linearity/Dynamic Range.....	23
	Assay Sensitivity.....	24

 Appendix.....	25
Bibliography.....	25
Cepheid Headquarters Locations	26
Technical Assistance.....	26
Table of Symbols	27
Revision History.....	28

Getting Started

Product Information

Proprietary Name

Xpert[®] BCR-ABL Ultra p190

Common or Usual Name

Xpert BCR-ABL Ultra p190

Product Description, Summary, and Principle of Procedure

Product Description

The Xpert BCR-ABL Ultra p190 assay, performed on the GeneXpert[®] Instrument Systems, is a real-time RT-PCR (reverse transcription-polymerase chain reaction) assay for the quantitative detection of the BCR-ABL1 p190 chromosomal translocation mRNA transcript (type e1a2) and the ABL1 endogenous control mRNA transcript in peripheral blood samples.

Summary and Explanation

The Philadelphia (Ph) chromosome is a shortened chromosome that results from the transposition of the 3' part of the ABL gene on chromosome 9 to the 5' part of the BCR gene on chromosome 22. The breakpoint on the ABL gene is fairly constant occurring at the 5' end of exon a2 whereas the breakpoint of the BCR gene is variable but is mainly clustered in 3 different regions (breakpoint cluster regions or bcr). Depending on the breakpoint on chromosome 22, different size segments are joined with the 3' sequences of the ABL gene. There are major (M-bcr), minor (m-bcr) and micro-breakpoints, each of which result in different size mRNA fusion transcripts.¹

The M-bcr results in the e13a2(b2a2) and e14a2(b3a2) fusion transcripts which are translated into BCR-ABL p210 proteins. The m-bcr results in the e1a2 fusion transcript and is translated into the BCR-ABL p190 protein. The e19a2 micro-breakpoint is translated into the BCR-ABL p230 protein.^{1,2}



The Ph chromosome is observed in greater than 95% of individuals with chronic myeloid leukemia (CML) and up to 20-30% of adults with acute lymphoblastic leukemia (ALL), 5% of children with ALL and in 1-2% of individuals with acute myeloid leukemia (AML).^{1,3}

In CML, the BCR-ABL p210 is present in greater than 95% of individuals and is also found in approximately 30% of Ph-positive (Ph+) ALL individuals. In the remaining individuals with Ph+ ALL and in rare cases of CML (1-3%), the BCR-ABL p190 is present.

In CML, the BCR-ABL p210 and p190 can co-exist. Both the p210 and p190 fusion proteins demonstrate increased tyrosine phosphokinase activity compared to the normal p145 c-abl protein.^{1,2,4,5}

In Ph+ ALL individuals, the p190 form is detected in approximately 80% of Ph+ childhood ALL and 20-40% of Ph+ adult ALL.¹

Principle of the Procedure

The Xpert BCR-ABL Ultra p190 assay is an automated assay for quantifying the amount of BCR-ABL1 p190 transcript as a ratio of BCR-ABL1 p190/ABL1. The assay is performed on Cepheid GeneXpert Instrument System, which automates and integrates sample purification, nucleic acid amplification, and target sequence detection in simple or complex samples using real-time RT-PCR and nested PCR assays. The system consists of an instrument, computer, and pre-loaded software for running assays and viewing the results. The system requires the use of single-use, disposable GeneXpert cartridges that hold the RT-PCR and nested PCR reagents and host the RT-PCR and nested PCR processes. For a full description of the system, refer to the relevant system operator manual.

The Xpert BCR-ABL Ultra p190 assay includes reagents to detect BCR-ABL1 p190 fusion genes resulting from a minor breakpoint, translocation e1a2, and the ABL1 transcript as an endogenous control in peripheral blood samples. The amount of BCR-ABL1 p190 transcript is quantified as the percent ratio of BCR-ABL1 p190/ABL1. There are two controls included in the Xpert BCR-ABL Ultra p190 assay – the Endogenous Control (ABL1) and a Probe Check Control (PCC). The ABL1 endogenous control normalizes the BCR-ABL1 p190 target and ensures that sufficient sample is used in the assay. The PCC verifies reagent rehydration, PCR tube filling, and that all reaction components, including probes and dyes, are present and functional in the cartridge.

Reagents, Instruments, and Materials

Reagents

Materials Provided

The Xpert BCR-ABL Ultra p190 kit (RBCRABLP190-10) contains sufficient reagents to process 10 test samples or quality control samples. The kit contains the following:

Xpert BCR-ABL Ultra Reagents	10 of each per kit
• Proteinase K (PK)	10 x 130 µL per vial
• Lysis Reagent (LY) (Guanidinium Chloride)	10 x 5.3 mL per vial
• Wash Reagent (1)	
◦ Ethanol	10 x 2.9 mL per ampoule
◦ Guanidinium thiocyanate	
Xpert BCR-ABL Ultra p190 Cartridges with Integrated Reaction Tubes	10 per kit
• Bead 1, 2, 3 and 4 (freeze-dried)	1 of each per cartridge
• Rinse Reagent	2.0 mL per cartridge
• Elution Reagent	2.5 mL per cartridge
CD	1 per kit
• Assay Definition File (ADF)	
• Instruction to import ADF into the GeneXpert software	
• Instructions for Use (Package Insert)	

Note Safety Data Sheets (SDS) are available at www.cepheid.com or www.cepheidinternational.com under



the **SUPPORT** tab.

Note The bovine serum albumin (BSA) in the beads within this product was produced and manufactured exclusively from bovine plasma sourced in the United States. No ruminant protein or other animal protein was fed to the animals; the animals passed ante- and post-mortem testing. During processing, there was no mixing of the material with other animal materials.

Note Certificates of Analysis and Lot Specifications Data Sheets are available through Cepheid Technical Support.

Materials Required but Not Provided

- GeneXpert system with touchscreen: GeneXpert instrument, touchscreen unit with built-in scanner, Cepheid OS software version 2.0 or higher, and operator manual.
- Printer: If a printer is required, contact Cepheid Technical Support to arrange for the purchase of a recommended printer.
- Vortex mixer
- Microcentrifuge (1,000 x g minimum)
- Pipettes and aerosol filter pipette tips
- 50 mL conical tubes
- Reagent grade absolute ethanol

Materials Available but Not Required

Xpert BCR-ABL Ultra p190 External Controls, INTROL[®] BCR-ABL1 p190 Control Panel, Catalog number C183, are quality controls from Maine Molecular Quality Controls, Inc.

Warnings and Precautions

General

- For Research Use Only. Not for use in diagnostic procedures. **RUO**
- Treat all biological samples, including used cartridges and reagents, as if capable of transmitting infectious agents. Because it is often impossible to know which might be infectious, all biological samples should be treated with standard precautions. Guidelines for sample handling are available from U.S. Centers for Disease Control and Prevention⁶ and Clinical and Laboratory Standards Institute.⁷
- Follow safety procedures set by your institution for working with chemicals and handling biological samples.
- The assay function has been established with blood collected in EDTA tubes only. The assay function has not been evaluated with other sample types.
- Reliable results are dependent on adequate sample collection, transport, storage and processing. Incorrect assay results may occur from improper sample collection, handling or storage, technical error, sample mix-up or because the target transcript in the sample is below the limit of detection of the assay. Careful compliance with the Package Insert instructions and the relevant system operator manual are necessary to avoid erroneous results.
- Performing the Xpert BCR-ABL Ultra p190 assay outside the recommended kit or sample storage temperature ranges and time may produce erroneous or invalid results.
- Biological samples, transfer devices, and used cartridges should be considered capable of transmitting infectious agents requiring standard precautions. Follow your institution's environmental waste procedures for proper disposal of used cartridges and unused reagents. These materials may exhibit characteristics of




chemical hazardous waste requiring specific national or regional disposal procedures. If national or regional regulations do not provide clear direction on proper disposal, biological samples and used cartridges should be disposed per WHO (World Health Organization) medical waste handling and disposal guidelines.⁸

Sample

- Maintain proper storage conditions during sample transport to ensure the integrity of the sample (see [Sample Collection, Transport and Storage](#), Sample Collection, Transport and Storage). Sample stability under shipping conditions other than those recommended has not been evaluated.
- Do not freeze whole blood samples.
- Proper sample collection, storage, and transport are essential for correct results.

Assay/Reagent


- Do not substitute Xpert BCR-ABL Ultra p190 reagents with other reagents.
- Do not open the Xpert BCR-ABL Ultra p190 cartridge lid except when adding sample and Wash Reagent.
- Do not use a cartridge that has been dropped after removing it from the packaging.



- Do not shake the cartridge. Shaking or dropping the cartridge after opening the cartridge lid may yield invalid results.
- Do not place the sample ID label on the cartridge lid or on the barcode label of the cartridge.



- Do not use a cartridge with a damaged barcode label.
- Hold the cartridge by the base. Do not touch the reaction tube at the rear of the cartridge, as this could cause damage that would interfere with light passing through it during the test. Do not use a cartridge that has a damaged reaction tube.



- It is recommended that the Xpert BCR-ABL Ultra p190 cartridges be at room temperature (20 °C – 30 °C) when used for testing.
- Each single-use Xpert BCR-ABL Ultra p190 cartridge is used to process one assay.
- Do not reuse processed cartridges.
- Do not reuse pipette tips.
- Do not use a cartridge if it appears wet or if the lid seal appears to have been broken.
- Do not use the Xpert BCR-ABL Ultra p190 cartridge if a reagent is added to the wrong opening.
- Do not open Xpert BCR-ABL Ultra p190 cartridges after the assay is completed.
- Dedicate a set of pipettes and reagents exclusively to sample preparation.
- Wear clean lab coats and gloves.
- Change gloves between the handling of each sample.
- In the event of a sample or control spill, wear gloves and absorb the spill with paper towels. Thoroughly clean and disinfect all laboratory work surfaces with a freshly prepared solution of 0.5% sodium hypochlorite in distilled or deionized water (dilute household bleach 1:10). Final active chlorine concentration should be 0.5%. After the work area is dry, follow by wiping the surface with 70% ethanol. For equipment, follow the manufacturer's recommendations for decontamination of equipment. Alternately, follow your institution's standard procedures for a contamination or spill event.

Chemical Hazards, Storage and Handling

Chemical Hazards^{9,10}

Note The information below applies to the Proteinase K, Lysis, Wash, and Rinse Reagents.

- UN GHS Hazard Pictogram:
- Signal Word: DANGER
- **UN GHS Hazard Statements**
 - Harmful if swallowed H302
 - Highly flammable liquid and vapor H225
 - Causes skin irritation H315
 - Causes serious eye irritation H319
 - May cause drowsiness or dizziness H336
 - Suspected of causing genetic defects H341
- **UN GHS Precautionary Statements**
 - **Prevention**



- Refer to Safety Data Sheet for special instructions before use.
- Do not handle until all safety precautions have been read and understood.
- Use personal protective equipment: gloves, eyewear, face shield and clothing.
- Use only in well-ventilated areas.
- Keep away from heat, sparks, open flames and/or hot surfaces.
- Avoid breathing mist, vapors, or spray.
- Wash hands thoroughly after handling.
- **Response**
 - In case of FIRE: Use appropriate media for extinction.
 - If INHALED: Remove victim to fresh air and keep at rest in a position comfortable for breathing.
 - Call a POISON CENTER or doctor/physician if victim feels unwell.
 - If SPILLED: Immediately remove contaminated clothing. If on skin or hair, rinse with water/shower.
 - If SKIN IRRITATION occurs: Get medical advice/attention.
 - If IN EYES: Remove contact lenses, if present. Rinse eyes thoroughly with water for several minutes. If eye irritation persists: Get medical advice/attention.
 - Specific treatment: see supplemental first aid measures in Safety Data Sheet.
 - If exposed or concerned: Get medical advice/attention.
- **Storage/Disposal**
 - Store under refrigerated conditions.
 - Keep containers tightly closed.
 - Dispose of content and/or container in accordance with local, regional, national, and/or international regulations.

Storage and Handling

- Store the Xpert BCR-ABL Ultra p190 kit contents at 2°C to 8°C until the expiration date provided on the label.
- Do not open the cartridge lid until you are ready to perform the assay.
- Do not use cartridges that have passed the expiration date.
- Do not use a cartridge that has leaked.
- The Wash Reagent is a clear, colorless liquid. Do not use the Wash Reagent if it has become cloudy or discolored.
- Twenty (20) minutes before starting the procedure, remove the blood sample, cartridge and sample preparation reagents from storage to allow them to come to room temperature (20°C to 30°C).

Specimen Collection, Testing, and Results

Specimen Collection



Specimen Collection, Transport and Storage

- Peripheral blood samples should be collected in EDTA tubes following your institution's guidelines. Plasma should not be separated from cells.
- Samples should be stored at 2°C to 8°C for no longer than 3 days (72 hours) prior to testing.
- Proper sample collection, storage, and transport are critical to the assay function. Sample stability under shipping and storage conditions other than those listed below have not been evaluated with the Xpert BCR-ABL Ultra p190 assay.

Procedure

Before You Start

Twenty (20) minutes before starting the procedure, remove the blood sample, sample preparation reagents, and cartridges from refrigerated storage to allow them to come to room temperature. Briefly spin down the Proteinase K (PK) in a microcentrifuge.

-  **Important** Remove the cartridge from the cardboard packaging before preparing the sample. (See **Preparing the Cartridge**).
-  **Important** Start the assay within 1 hour of adding the Sample Reagent-treated sample to the cartridge.

Preparing the Sample

Preparing a Sample with Unknown White Blood Cell (WBC) Count or with No Greater than 30 Million WBC/mL

1. To the bottom of a new 50 mL conical tube, add 100 μ L of PK (Proteinase K).
2. Ensure blood sample is well-mixed by inverting the blood collection tube 8 times immediately before pipetting. See manufacturer's instructions for the EDTA blood collection tube.



3. To the tube already containing Proteinase K, add 4 mL of blood sample.
4. Mix the sample with a vortex mixer at maximum setting continuously for 3 seconds.
5. Incubate at room temperature for 1 minute.
6. To the same tube, add 2.5 mL of Lysis Reagent (LY).
Note Retain the remaining lysis reagent to use again in Step 13.
7. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds.
8. Incubate at room temperature for 5 minutes.
9. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds.
10. Incubate at room temperature for 5 minutes.
11. Mix the sample by tapping the bottom of the tube 10 times.
12. Transfer 1 mL of the prepared lysate into a new 50 mL conical tube.
Note Remaining lysate can be stored at 2–8 °C for up to 48 hours or stored at -20 °C or lower for up to 3 months.
13. To the new conical tube containing lysate, add 1.5 mL of retained Lysis Reagent (LY) from Step 6.
14. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds.
15. Incubate at room temperature for 10 minutes.
16. To the same conical tube, add 2 mL of reagent grade absolute ethanol (provided by user).
17. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds. Set aside.
18. Discard any remaining PK or LY reagents.

Preparing a Sample with WBC Count Greater than 30 Million cells/mL

1. To the bottom of a new 50 mL conical tube, add 100 µL of PK (Proteinase K).
2. Ensure blood sample is well-mixed by inverting the blood collection tube 8 times immediately before pipetting. See manufacturer's instructions for the EDTA blood collection tube.
3. To the tube already containing Proteinase K, add 50 µL of blood sample.
4. Mix the sample with a vortex mixer at maximum setting continuously for 3 seconds.
5. Incubate at room temperature for 1 minute.
6. To the same tube, add 2.5 mL of Lysis Reagent (LY).
7. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds.
8. Incubate at room temperature for 5 minutes.
9. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds.
10. Incubate at room temperature for 5 minutes.
11. To the same conical tube, add 2 mL of reagent grade absolute ethanol (provided by user).
12. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds. Set aside.
13. Discard any remaining PK or LY reagents.

Preparing the Cartridge

To add the sample to the Xpert BCR-ABL Ultra p190 cartridge:



1. Remove the cartridge from the cardboard packaging.
2. Inspect the cartridge for damage. Do not use if damaged.
3. Lift the cartridge lid and transfer the entire contents of the Wash Reagent (1) ampoule to the Wash Reagent Chamber (small opening). See [Figure 1](#).
4. Pipette the entire contents of the prepared sample into the Sample Chamber (large opening). See [Figure 1](#).




Figure 1 Xpert BCR-ABL Ultra p190 Cartridge (Top View)

5. Close the cartridge lid. Ensure the lid snaps firmly into place. Initiate assay (see Starting the Test).


Starting the Test: GeneXpert System with Touchscreen

i Important Before you start the test, make sure that:

- The system is running the correct Cepheid OS software version shown in section - **Materials Required but Not Provided**.
- The correct assay definition file is imported into the software.


Note The default workflow is shown. Your system administrator may alter the workflow.

1. Turn on GeneXpert system with touchscreen.
2. Log on to system software using your username and password.
3. On the Modules tab, touch **Start Test**.
4. Follow onscreen prompts to create new test and enter patient and sample information.
5. Scan or manually input the cartridge serial number. If scanning, hold the cartridge about 1-3 inches (3-7 cm) away from the scanner. The scanner projects a green crosshair, which you center on the barcode. Scanning is complete when you hear an audible beep. Touch **Continue**.






6. Select the desired test and touch **Continue**.
7. Watch the cartridge preparation video, if needed.
8. On the Confirm screen, review all data and touch **Confirm**.
9. Open the module door under flashing green light and insert the cartridge.



10. Close cartridge module door completely by pressing until it latches. The test starts.



11. When the test completes, the **Results Summary** screen appears. Open the module door and remove cartridge.
12. Dispose of used cartridge in appropriate waste container according to your institution's standard practices.

Viewing Results: GeneXpert System with Touchscreen

The GeneXpert system with touchscreen results screen will automatically interpret test results for you and clearly show them in the **View Results** window.

1. Tap **Results**.
2. Tap the test to be viewed in the Results screen.
3. Click **OK**.
4. To generate a PDF report file, touch **View Report**. More detailed instructions for viewing and uploading results are available in your system operator manual.

Quality Control

Each cartridge includes an ABL 1 Endogenous Control and a Probe Check Control (PCC).



ABL 1 Endogenous Control — The ABL 1 Endogenous Control verifies that sufficient sample is used with the assay. Additionally, this control detects sample-associated inhibition of the real-time PCR assay. The ABL 1 passes if it meets the assigned acceptance criteria.

Probe Check Control (PCC) — Before the start of the PCR reaction, the GeneXpert system measures the fluorescence signal from the probes to monitor bead rehydration, reaction tube filling, and if all reaction components are functional in the cartridge. The PCC passes if it meets the assigned acceptance criteria.

Results

Xpert BCR-ABL Ultra p190 quantitative outputs are provided as a percent ratio of BCR-ABL1 p190/ABL1. Examples of possible results and interpretations are presented in [Table 1](#).

Table 1. Xpert BCR-ABL Ultra p190 Possible Results and Interpretation

Probe Check*	ABL Ct*	e1a2 Ct*	Xpert BCR-ABL Ultra p190 Test Result	Notes
PASS	PASS	POS	BCR-ABL p190 DETECTED [#.###%]	Calculated% ratio value is reported.
			BCR-ABL p190 DETECTED [Above upper LoQ]	Calculated% value is above the limit of quantitation and is not reported.
			BCR-ABL p190 DETECTED [Below LoD;<0.0065%]	Calculated% ratio is below the limit of detection and is not reported.
		NEG	BCR-ABL p190 NOT DETECTED [Sufficient ABL transcript]	e1a2 Ct is zero or above the acceptance threshold.
		INVALID	INVALID [Too high BCR-ABL p190 transcript]	e1a2 Ct is below the acceptance threshold.
		FAIL	POS, NEG, or INVALID	INVALID [No ABL transcript]
	INVALID [Insufficient ABL transcript]			ABL Ct is above the acceptance threshold.
	INVALID [Too high ABL transcript]			ABL Ct is below the acceptance threshold.
	INVALID		INVALID [Too high BCR-ABL p190 and ABL transcripts]	Both e1a2 and ABL Ct values are below the acceptance thresholds.
	FAIL	PASS or FAIL	POS, NEG or INVALID	ERROR

* See the Analyte Results tab in the GeneXpert System Software for details

Quantitative Results

Xpert BCR-ABL Ultra p190 quantitative outputs are provided as a percent ratio of BCR-ABL1 p190/ABL1. Kits are assigned lot-specific Efficiency (EΔCt) and Scaling Factor (SF) values that tie the quantitation of BCR-ABL1 p190 (e1a2) and ABL1 transcripts to copy numbers of synthetic BCR-ABL p190 and ABL1 RNA in vitro transcribed RNA (IVT-RNA) primary standards.



Table 2. Examples of Xpert BCR-ABL Ultra p190 Assay Results

Test	BCR-ABL p190		ABL		Xpert BCR-ABL Ultra p190 Assay Results	Notes
	Ct	Result	Ct	Result		
1	7.1	INVALID	7.3	FAIL	INVALID [Too high BCR-ABL p190 and ABL transcripts]	Calculated % value: 203.17%
2	8.1	INVALID	7.9	FAIL	INVALID [Too high ABL transcript]	Calculated % value: 152.46%
3	7.9	INVALID	8.1	PASS	INVALID [Too high BCR-ABL p190 transcript]	Calculated % value: 203.17%
4	25.0	INVALID	18.2	FAIL	INVALID [Insufficient ABL transcript]	NA
5	0	INVALID	0	FAIL	INVALID [No ABL transcript]	NA
6	12.3	POS	11.2	PASS	BCR-ABL p190 DETECTED [Above upper LoQ]	Calculated % value: 79.91%
7	22.5	POS	11.8	PASS	BCR-ABL p190 DETECTED [0.081%]	Calculated % value: 0.081%
8	27.8	POS	11.6	PASS	BCR-ABL p190 DETECTED [Below LoD; <0.0065%]	Calculated % value: 0.0016%
9	0	NEG	12.1	PASS	BCR-ABL p190 NOT DETECTED [Sufficient ABL transcript]	0%
10	0	NO RESULT	0	NO RESULT	ERROR	For example, Error 5017 [ABL] probe check failed

Note

GeneXpert systems calculate results automatically based upon *cycle threshold* (Ct) values generated by the test, and lot-specific parameters assigned during manufacturing. The software applies the following algorithm, wherein the Δ Ct (Delta Ct) value is obtained from ABL Ct minus BCR-ABL p190 Ct, and Efficiency (E) and Scaling Factor (SF) are lot specific values:

$$\text{Percent ratio} = \text{Efficiency}^{(\Delta\text{Ct})} \times \text{Scaling Factor} \times 100$$

Note Efficiency and Scaling Factor values calibrate the quantitation of BCR-ABL1 p190 (e1a2) and ABL1 transcripts to copy numbers of synthetic BCR-ABL p190 and ABL1 RNA *in vitro* transcribed RNA (IVT-RNA) primary standards. Efficiency and Scaling Factor values are embedded within each cartridge barcode. Lot Specifications Data Sheets are available through Cepheid Technical Support.

BCR-ABL p190 DETECTED [#.##%]

For a “BCR-ABL p190 DETECTED [#.##%]” result, BCR-ABL p190 is detectable with BCR-ABL p190 Ct greater than or equal to “8” and less than or equal to the cut-off of “32” and ABL Ct greater than or equal to “8” and less than or equal to “18”.



Lot-specific $E_{\Delta Ct} = 2.05$; $SF = 1.76$

Example: Assay's ABL Ct = 11.4; BCR-ABL p190 Ct = 15.6; $\Delta Ct = -4.2$

$$\% = 2.05^{(-4.2)} \times 100 \times 1.76 = 8.63\%$$

Result: BCR-ABL p190 DETECTED [8.63%].

BCR-ABL p190 DETECTED [Above upper LoQ]

BCR-ABL p190 has been detected at a level > 25%.

For a “BCR-ABL p190 DETECTED [Above upper LoQ]” result, BCR-ABL p190 is detectable with BCR-ABL p190 Ct greater than or equal to “8” and less than or equal to the cut-off of “32” and ABL Ct greater than or equal to “8” and less than or equal to “18”.

Lot-specific $E_{\Delta Ct} = 2.05$; $SF = 1.76$

Example: Assay's ABL Ct = 17.2; BCR-ABL p190 Ct = 18.7; $\Delta Ct = -1.6$

$$\% = 2.05^{(-1.6)} \times 100 \times 1.76 = 56.6\% \text{ is greater than the defined assay upper LoQ at } 25\%$$

Result: BCR-ABL p190 DETECTED [Above upper LoQ].

BCR-ABL p190 DETECTED [Below LoD; <0.0065%]

BCR-ABL p190 has been detected at a level < 0.0065%.

For a “BCR-ABL p190 DETECTED [Below LoD; <0.0065%]” result, BCR-ABL p190 is detectable with BCR-ABL p190 Ct greater than or equal to “8” and less than or equal to the cut-off of “32” and ABL Ct greater than or equal to “8” and less than or equal to “18”.

Lot-specific $E_{\Delta Ct} = 2.05$; $SF = 1.76$

Example: Assay's ABL Ct = 10.1; BCR-ABL p190 Ct = 24.8; $\Delta Ct = -14.8$

$$\% = 2.05^{(-14.8)} \times 100 \times 1.76 = 0.0044\% \text{ is less than the defined assay LoD at } 0.0065\%$$

Result: BCR-ABL p190 DETECTED [Below LoD; <0.0065%].

BCR-ABL p190 NOT DETECTED [Sufficient ABL transcript]

BCR-ABL p190 was not detected with BCR-ABL p190 Ct equal to “0” or greater than the cut-off of “32” and ABL Ct greater than “8” and less than or equal to “18”.

When BCR-ABL p190 is undetectable with BCR-ABL p190 Ct equal to “0” or greater than the cut-off of “32”, the GeneXpert software first looks for the ABL Ct to confirm if the ABL Ct is greater than or equal to “8” and less than or equal to “18” to ensure having “Sufficient ABL transcript”. See [Table 2](#).

Example: Assay's BCR-ABL p190 Ct = 0; ABL Ct = 11.6 is less than “18”.

Result: BCR-ABL p190 NOT DETECTED [Sufficient ABL transcript].

INVALID [No ABL transcript]

BCR-ABL p190 was not detected with ABL Ct equal to “0”.

When BCR-ABL p190 is either detected or not detected, the GeneXpert software first looks for the ABL Ct to confirm if the ABL Ct is less than or equal to “18” to ensure having “Sufficient ABL transcript”. Refer to [Troubleshooting Guide](#).

Example: Assay's BCR-ABL p190 Ct = 0; ABL Ct = 0.

Result: INVALID [No ABL transcript].



INVALID [Insufficient ABL transcript]

BCR-ABL p190 was not detected with ABL Ct greater than "18".

When BCR-ABL p190 is either detected or not detected, the GeneXpert software first looks for the ABL Ct to confirm if the ABL Ct is less than or equal to "18" to ensure having "Sufficient ABL transcript". Refer to [Troubleshooting Guide](#).

Example: Assay's BCR-ABL p190 Ct = 31.2; ABL Ct = 28 is greater than "18".

Result: INVALID [Insufficient ABL transcript].

INVALID [Too high BCR-ABL p190 and ABL transcripts]

BCR-ABL p190 was detected with both BCR-ABL p190 and ABL Cts less than "8".

When BCR-ABL p190 is either detected or not detected, the GeneXpert software first looks for the ABL Ct to confirm if the ABL Ct is less than or equal to "18" to ensure having "Sufficient ABL transcript". Refer to [Troubleshooting Guide](#).

Example: Assay's BCR-ABL p190 Ct = 7.9; ABL Ct = 7.6 is less than "8".

Result: INVALID [Too high BCR-ABL p190 and ABL transcripts].

ERROR

Refer to [Troubleshooting Guide](#).

Troubleshooting Guide

Table 3. Troubleshooting Guide

Assay Result	Possible Causes	Suggestions
INVALID	Type 1: Endogenous control ABL failure: <ul style="list-style-type: none"> Poor sample quality RT-PCR inhibition If ABL Ct > 18, and/or endpoint <200 	<ul style="list-style-type: none"> Check the sample quality (e.g., exceeded sample storage requirement including time and temperature). Repeat the assay with original sample (if available) or from retained lysate and a new cartridge following the procedure as described in , Retest Procedure for ERROR or INVALID (Type 1).
	Type 2: BCR-ABL p190 transcript level cannot be determined due to sample containing excess BCR-ABL p190 and/or ABL transcripts (Ct < 8)	Repeat the assay with original sample (if available) or from retained lysate and a new cartridge following the procedure as described in , Retest Procedure for ERROR (Code 2008) or INVALID (Type 2).
ERROR (Code 2008)	Pressure exceeding limit (error message 2008)	<ul style="list-style-type: none"> Check the sample quality Check for grossly elevated WBC count Repeat the assay with original sample (if available) or from retained lysate and a new cartridge following the procedure as described in , Retest Procedure for ERROR (Code 2008) or INVALID (Type 2).



Assay Result	Possible Causes	Suggestions
ERROR (Code 5006, 5007, 5008, and 5009 ^a)	Probe check failure	Repeat the assay with original sample (if available) or from retained lysate and with a new cartridge following the procedure as described in , Retest Procedure for ERROR or INVALID (Type 1).
NO RESULT	Data collection failure. For example, the operator stopped an assay that was in progress or a power failure occurred.	Repeat the assay with original sample (if available) or from retained lysate and with a new cartridge following the procedure as described in , Retest Procedure for ERROR or INVALID (Type 1).

a. This is not an exhaustive list of ERROR codes.

Retests

Retest Procedure for ERROR or INVALID (Type 1)

Retest samples with **ERROR** or **INVALID** results due to the ABL cycle threshold (Ct) exceeding the maximum valid Ct cut-off (Ct >18) or the endpoint is below the threshold setting (< 200). Also refer to [Troubleshooting Guide](#).

1. Measure blood sample volume:

- If *sufficient* blood sample volume is available, re-test from original blood sample collection tube following the procedure in [Preparing the Sample](#).

-OR-

- If blood sample volume is insufficient, re-test can be performed with the retained lysate from [Preparing the Sample](#), Step 12.
 - a. If retained lysate from [Preparing the Sample](#), Step 12 is stored frozen, thaw to room temperature before use.
 - b. Ensure lysate is well-mixed by mixing the sample with a vortex mixer at maximum setting continuously for 10 seconds and set it aside for 3 minutes for bubbles to settle. Go to Step 2.

2. Transfer 1 mL of the prepared lysate into a new 50 mL conical tube.

3. To the new conical tube containing lysate, add 1.5 mL of Lysis Reagent (LY).

4. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds.

5. Incubate at room temperature for 10 minutes.

6. To the same conical tube, add 2 mL of reagent grade absolute ethanol (not provided).

7. Mix the sample with a vortex mixer at maximum setting continuously for 10 seconds.

8. Open the cartridge by lifting the cartridge lid and transfer the entire contents of the Wash Reagent (1) ampoule to the Wash Reagent chamber (with small opening). See [Figure 1](#).

9. Pipette the entire contents of the prepared sample into the Sample Chamber (large opening). See [Figure 1](#).

10. Close cartridge lid. Initiate assay (see Starting the Test).

Retest Procedure for ERROR (Code 2008) or INVALID (Type 2)

Retest samples with BCR-ABL p190 and/or ABL transcript levels below the valid minimum Ct cut-off (Ct < 8) and/or when pressure limit is exceeded. Also refer to [Troubleshooting Guide](#).

1. To the bottom of a new 50 mL conical tube, add 100µL of PK (Proteinase K).



2. Measure blood sample volume:

- If *sufficient* blood sample volume is available, re-test from original blood sample collection tube. Ensure blood sample is well-mixed by inverting the blood collection tube 8 times immediately before pipetting. Go to Step 4.

-OR-

- If blood sample volume is *insufficient*, re-test can be performed from the retained lysate from [Preparing the Sample](#) Step 12.
 - a** If retained lysate from [Preparing the Sample](#) Step 12 is stored frozen, thaw to room temperature before use. If refrigerated lysate is used, allow to come to equilibrate to room temperature before use.
 - b** Ensure lysate is well-mixed by mixing the sample with a vortex mixer at maximum setting continuously for 10 seconds and set it aside for 3 minutes for bubbles to settle. Go to Step 3.
3. To the tube already containing Proteinase K, add 50 µL of blood sample, if available, or 80 µL of left-over lysate from [Preparing the Sample](#).
 4. Mix the sample with a vortex mixer at maximum setting continuously for 3 seconds.
 5. Incubate at room temperature for 1 minute.
 6. Follow the Steps 3-13 in to make the final lysate.
 7. Open the cartridge by lifting the cartridge lid and transfer the entire contents of the Wash Reagent (1) ampoule to the Wash Reagent chamber (with small opening). See [Figure 1](#).
 8. Pipette the entire contents of the prepared sample into the Sample Chamber (large opening). See [Figure 1](#).
 9. Close cartridge lid. Initiate assay (see Starting the Test).

Limitations

Limitations of the Assay

- For Research Use Only. Not for use in diagnostic procedures.
- The assay is not intended to be used with external calibrators.
- The performance of the Xpert BCR-ABL Ultra p190 RUO assay was evaluated using the procedures provided in these Instructions for Use only. Modifications to these procedures may alter the performance of the assay.
- This product was designed for use with blood collected in EDTA tubes only.
- Do not use heparin as the anticoagulant because it can inhibit the PCR reaction.
- Sodium citrate, buffy-coat and bone marrow sample types have not been validated.
- Xpert BCR-ABL Ultra p190 RUO assay is designed to detect the p190 BCR-ABL fusion transcript e1a2. The assay does not detect major or micro breakpoints, microdeletions, or mutations.
- The ability to detect other fusion transcripts has not been evaluated beyond those described in these instructions for use. The Xpert BCR-ABL Ultra p190 assay is not intended to detect the e13a2/b2a2 and e14a2/b3a2 (p210), e19a2 (p230) or other minor translocations.
- Samples with high white blood cell counts (greater than 30 million cells/mL) tested per may yield inaccurate results, **INVALID** results, or aborted runs due to pressure build-up within the cartridge. See [Table 3](#) for additional information.
- Samples with very low levels of ABL transcript or with white blood cells lower than 150,000 cells/mL may be reported as **INVALID** (Type 1). A non-determinate result does not preclude the presence of very low levels of leukemic cells in the sample.



- Mutations or polymorphisms in primer or probe binding regions may affect detection of new or unknown variants and may result in a false negative result.
- Some samples with very low levels of BCR-ABL1 transcript (i.e., below LoD 0.0065%) may be reported as **BCR-ABL p190 NOT DETECTED [Sufficient ABL transcript]**. Hence, an undetected result does not preclude the presence of low levels of leukemic cells in the sample.

! Specific Performance Characteristics

Analytical Data

The data was collected from internal studies only.

Assay Linearity/Dynamic Range

The linearity of the assay was evaluated by serially diluting total RNA isolated from SUP-B15 (e1a2) cell line in a background lysate prepared from normal EDTA peripheral blood samples to target 9 different levels ranging between 25% to 0.001% BCR-ABL p190/ABL and testing each level, including the negative, on two assay kit lots for a total of 80 replicates.

Linear regression analyses were performed for first, second and third order polynomials. The results were considered linear if the polynomial regression coefficients were insignificant (p -values > 0.05). The combined linear regression curve for e1a2 transcript from two kit lots is shown in [Figure 2](#).




Figure 2 Linear Regression Curves for p190 Transcript (e1a2)



The estimated regression intercept, slope, and R^2 value from the linear model are shown in [Table 4](#).

Table 4. Regression Coefficients from Linear Model

Intercept	Slope	R^2
-0.0561	0.9248	0.9811

Assay Sensitivity

The analytical limit of blank (LoB) was established using 40 normal EDTA peripheral blood samples. Two out of 130 total replicates tested reported as BCR-ABL p190 **DETECTED**. The overall LoB was determined to be 0.00031%.

The analytical limit of detection (LoD) was estimated for the p190 transcript using dilutions of SUP-B15 cell line in negative EDTA blood background. The estimated analytical LoD is the lowest analyte concentration (reported as % BCR-ABL p190/ ABL transcript) that can be reliably distinguished from the LoB > 95% of the time. The estimated analytical LoD of the assay is 0.0065% BCR-ABL p190/ABL.

The limit of quantitation (LoQ) was estimated with the data obtained from the LoD study. The mean for the % values was calculated for replicates at levels equal to the LoD, 0.0065%, or greater with positivity greater or equal to 95%. The LoQ of the assay is constrained by the LoD of the assay; therefore, the LoQ was determined to be equal to the LoD, 0.0065%.

? Appendix

Bibliography

1. Faderl S. et al. Clinical Significance of Cytogenic Abnormalities in Adult Acute Lymphoblastic Leukemia. *Blood*. 1998; 91 (11): 3995-4019.
2. Dushyant, V. et al. Chronic myeloid leukemia (CML) with p190 BCR-ABL: analysis and characteristics, outcomes and prognostic significance. *Blood*. 2009; 114: 2232-2235.
3. Kurzrock, R. et al. The Molecular Genetics of Philadelphia Chromosome-Positive Leukemias. Mechanisms of Disease. *NEJM* 1998; 990-998.
4. Chan, L.C. et al. A novel abl protein expressed in Philadelphia chromosome positive acute lymphoblastic leukemia. *Nature*. 1987; 325: 635-637.
5. Deininger M.W. The molecular biology of chronic myeloid leukemia. *Blood*. 2000; 96: 3343-3356.
6. Centers for Disease Control and Prevention. Biosafety in Microbiological and Biomedical laboratories (refer to latest edition).
7. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline. Document M29 (refer to latest edition).
8. <https://>
9. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.
10. Occupational Safety and Health Standards, Hazard Communication, Toxic and Hazard Substances (March 26, 2012) (29 C.F.R., pt. 1910, subpt. Z).



Cepheid Headquarters Locations

Corporate Headquarters

Cepheid
904 Caribbean Drive
Sunnyvale, CA 94089
USA

Telephone: + 1 408 541 4191
Fax: + 1 408 541 4192
www.cepheid.com

European Headquarters

Cepheid Europe SAS
Vira Solelh
81470 Maurens-Scopont
France

Telephone: + 33 563 825 300
Fax: + 33 563 825 301
www.cepheidinternational.com

Technical Assistance

Before contacting Cepheid Technical Support, collect the following information:

- Product name
- Lot number
- Serial number of the instrument
- Error messages (if any)
- Software version and, if applicable, Computer Service Tag number

United States Technical Support

Telephone: + 1 888 838 3222
Email: techsupport@cepheid.com

Contact information for all Cepheid Technical Support offices is available on our website: www.cepheid.com/en/support/contact-us.

Table of Symbols

Symbol	Meaning
	Catalog number
	Research Use Only
	Batch code
	Do not reuse
	Consult instructions for use
	Manufacturer
	Country of manufacture
	Contains sufficient for <i>n</i> tests
	Control
	Expiration date
	Temperature limitation
	Biological risks
	Caution
	Warning



Cepheid
904 Caribbean Drive
Sunnyvale, CA 94089
USA

Phone: + 1 408 541 4191

Fax: + 1 408 541 4192

www.cepheid.com



Revision History

Description of Changes: 303-0947, Rev. A

Purpose: Initial release

